

Proteomic sequencing and resurrection of a monoclonal antibody

Natalie Castellana¹, Kexin Huang², Jingxing Li², and Hua Tu²

¹ Digital Proteomics LLC, San Diego, CA, USA

² Curia, Inc., Belmont, CA, USA



Motivation

While hybridoma sequencing has become a routine part of drug development workflows, there are few options for antibody primary sequence determination when the source hybridoma is unavailable. Mass spectrometry offers a high-throughput method for interrogating the proteome. Digital Proteomics developed Valens, which combines an antibody-specific mass spectrometry protocol with sophisticated informatics to address the challenge of proteomic sequencing of antibodies.

FLAG[®] tag is a polypeptide tag that can be conjugated to a protein, enabling quantitative investigation of a protein(s) of interest. The ANTI-FLAG[®] M2 antibody targets FLAG tag and is a crucial component of the proteomic assays. A preliminary protein sequence of an ANTI-FLAG tag antibody was recently published for the purposes of elucidating the crystal structure of the antibody in complex with trans-membrane proteins (Roosild *et al*, 2006).

Protein sequencing was performed using Edman degradation; however, this method only revealed the first 40 amino acids and several gaps persisted, including in the CDR3. Recently, Digital Proteomics developed a monoclonal antibody sequencing technology, Valens, which employs mass spectrometry to determine the full antibody sequence. In collaboration with Curia, Digital Proteomics determined the antibody sequence for the ANTI-FLAG M2 antibody.

Protein digestion

A single MS run (~4500 spectra) of a single enzyme only covers about 50% of the antibody sequence. Using our recommended five enzymes (trypsin, chymotrypsin, pepsin, elastase, and AspN) routinely achieves over 99% coverage.

QIQLQQSGAELASPGASVKLSCKASGYTFTDHIMNWWKKRPGQGLEWIGRIYPVKKETNY		
QIQLQQSGAELASPGASVK		
	ASGYTFTDHIMNWK	
		RPGQGLEWIGR
LQQSGAEL		
	SCKASGY	
		VKKRPGQGLEW
QIQLQQSGAELASPGASVKLSCKASGYTFT		
		DHIMNWWKKRPGQGLEWIGRIYPVKKETNY
Trypsin	Chymotrypsin	AspN

Final sequences

Heavy Chain	
Valens	EVQLQQSAELARPGASVKMSCKASGYFTXIXIHWVKQRPGQGLEWIGYIPXXGXXXYNQNFKDETTLTADPSSS
Published	EVQLQQSGGEL?KPGASVKMSCK?SGYTFT?Y?IHW?KQ?-G?GLEWIGYI?P??G??-YN??FKGK?TL??DKSSS
Valens	TAYMELNSLTSEDSAVYYCARXXXXXDYWGQGLTLVSSAKTTPPSVYPLAPGSAAQTNSMVTLGCLVKGYFPEPV
Published	TAYM????LTSEDSAVY-C?R????G?DYWGQGLTLVSSAKTTPPSVYPLAPGSAAQTNSMVTLGCLVKGYFPEPV
Valens	TVTWNSGSLSSGVHTFPAVLQSDLYTLSSSVTVPSSTPWSETVTCNVAHPASSTKVDKKIV . . .
Published	TLTWNSGSLSSGVHTFPAVLQSDLYTLSSSVTVPSSTPWSETVTCNVAHPASSTKVDKKIV
Light Chain	
Valens	DVLMTIQIPLSLPVSLGDQASISCRSSQXIVHXNGNTYLEWYLLKPGQSPKLLIYKVNRFSGVPDRFSGSGSDTFT
Published	DVLMTIQ?PLSLPVSLGDQASISCRSSQ?IVH?NGNTYLEWYLQKPGQSP?LLIYKV?NRFSGVPDRFSGSGSDTFT
Valens	LKISRVEAEDLGVYYCFQGXHPYTFGGGKLEIKRADAAPTIVSIFPPSSEQLTSGGASVVCFLNFPKDIINVKWK
Published	LKISRVEAEDLGVYYCFQGH?PYTFGGGKLEIKRADAAPTIVSIFPPSSEQLTSGGASVVCFLNFPKDIINVKWK
Valens	IDGSERQNGVLSWTDQSKDSTYSMSSTLTLTKDEYERHNSYTCEATHKTSTSPIVKSFNR . . .
Published	IDGSERQNGVLSWTDQSKDSTYSMSSTLTLTKDEYERHNSYTCEATHKTSTSPIVKSFNR

34 unknown residues in the the published sequence ([periwinkle](#))


14 amino acid differences ([plum](#))

Some amino acids are replaced by 'X' for proprietary reasons.

Making recombinant HM2 antibody

A recombinant antibody with human constant regions and the amino acid sequences identified by Valens as the variable regions were made. Plasmid DNA encoding the full-length antibody was transfected to CHO cells to produce the HM2 antibody. HM2 antibody can be purified very


efficiently from the conditioned media with Protein A beads. Purified HM2 antibody was tested for binding to multiple proteins in comparison to Sigma ANTI-FLAG® M2 antibody, and the results showed that they have similar binding profile and affinity (Kd).




	Kd against FLAG tagged protein (Octet® system)	Kd against His-tagged protein (Octet® system)	Kd against 3x FLAG tagged protein (Biacore® system)
HM2	26 ± 3 nM	No binding	< 1 nM
M2 (reference)	25 ± 4 nM	No binding	< 1 nM

Using HM2 for purification

Input




Elute



HM2 antibody was chemically conjugated to sepharose beads. The resulting beads were used to purify from a mixture containing a FLAG tagged protein. As shown in the graphs above, the HM2 beads can specifically and efficiently purify FLAG

tagged proteins, with a recovery rate greater than 90%. The target protein binding capacity is 13.5 mg per mL of HM2 resin, which is at least 10x better than Sigma M2 resin.

Proteomic sequencing using Valens





ABOUT CURIA

Curia is a Contract Development and Manufacturing Organization with over 30 years of experience, an integrated network of 29 global sites and over 3,500 employees partnering with customers to make treatments broadly accessible to patients. Our biologics and small molecule offering spans discovery through commercialization, with integrated regulatory and analytical capabilities. Our scientific and process experts and state-of-the-art facilities deliver best-in-class experience across drug substance and drug product manufacturing. From curiosity to cure, we deliver every step to accelerate and sustain life-changing therapeutics. ***Learn more at curiaglobal.com***

Solutions developed by Curia and Digital Proteomics

CONTACT US

www.curiaglobal.com

