



PLASMIDS FOR PROTEIN OR ANTIBODY PRODUCTION

Clear Form

SENDER AND PROJECT INFORMATION:

Form with fields: Date, Reference, Name, Institution, Phone #, Email. Includes a note: Quote ID or Catalog #

Other information:

Hazards: [] Not Hazardous [] Biohazard (Attach MSDS if available)

Options (select one):

- Options 1, 2, and 3 regarding plasmid DNA scale up, sequence confirmation, and transfection-ready DNA requirements.

Main data table with columns: Construct Name/ ID, # vials, Concentration, Vol., Bacterial Resistance, Solvent, Endotoxin Level, Species, Isotype, Vector.

Ship samples to:

Curia Bio, Inc.
Attn. Receiving
21021 Corsair Blvd
Hayward, CA 94545
650-288-4891

The following QC criteria is typically required for Curia productions, and is strongly recommended for client-provided DNA ready for transfection:

- QC criteria list: DNA amount(s), Sequence confirmation, A260/A280, A260/230, Concentration, ccc/OC ratio, and filtering requirements.

* For options 1 or 2, please attach sequence or sequencing primers